

# Freeze-Pak™ STS Secondary Shipping Containers

The Freeze-Pak™ STS single-use secondary shippers are engineered to support frozen storage and transport needs intended for temperatures as low as -80°C. The STS shippers have been designed to complement the single-use Freeze-Pak™ STS bio-containers. They provide a simple solution while delivering the flexibility and durability required for supporting frozen storage and transport.


- Single-use
- Durable
- Frozen storage
- Frozen transport
- Compatible with Freeze-Pak™ STS single-use bags



## STS Shipper Specifications


### STS-XS

- Dimensions (fully assembled)
- 19.6" L x 11.8" W x 4.75" D  
49.8 cm x 29.8 cm x 12.1 cm
- Materials
  - Case: polypropylene plastic corrugate; blue
  - Insert: folding polypropylene plastic corrugate with polyurethane film



### STS-XL

- Dimensions (fully assembled)
- 26.8" L x 15" W x 4.9" D  
68.1 cm x 38.1 cm x 12.4 cm
- Materials
  - Case: polypropylene plastic corrugate; blue
  - Insert: folding plastic corrugate with polyurethane film
  - Sliding tray/pad: 5mm polypropylene plastic corrugate
  - Foam pad: Formed Instapak



## Ordering Information

Catalog Number	Description	Units per Case
STS-XS	Single-use, fully assembled secondary container intended for Freeze-Pak™ STS bags up to 1L	5
STS-XL	Single-use, fully assembled secondary container intended for Freeze-Pak™ STS bags up to 5L	5

Freeze-Pak™ is a trademark of charter Medical, Ltd. All data and statements concerning these products may be considered as being indicative of representative properties and characteristics obtainable. We make no warranty, express or implied, concerning actual use or results because of industry specific influences.



# CharterMedical

3948-A Westpoint Blvd. | Winston-Salem, NC 27103, USA | Toll Free: 866.458.3116 | Fax: 336.714.4241 | [www.chartermedical.com](http://www.chartermedical.com)

FP-ST5-SSC-2015-REV1