



Cryopreservation Media

INCELL Cryopreservation Media is manufactured using the highest quality components. Each formulation is chemically defined and animal-component free. The two primary product families are EZ-CPZ™ and StorZ™. EZ-CPZ™ and EZ-CPZ-ND™ products are designed for use with cell suspensions, whereas StorZ™ products are intended for use with cell pellets or tissues. Each of these media are available with and without DMSO. Solutions containing DMSO are light sensitive. All products must be stored at 2°C to 8°C and should not be frozen prior to use.

EZ-CPZ™

- » Ideal for use with suspended cells or for tissue storage
- » Contains 10% DMSO
- » 2X solution

EZ-CPZ-ND™

- » Ideal for use with suspended cells or for tissue storage
- » Without DMSO
- » 2X solution

StorZ™

- » Ideal for use with cell pellets or tissues
- » Contains 7% DMSO
- » 1X solution

StorZ-ND™

- » Ideal for use with cell pellet or tissues
- » Without DMSO
- » 1X solution

Ordering Information

Catalog #	Description	Case Qty.
EZCZ-100	100 mL EZ-CPZ™	1
EZCZ-500	500 mL EZ-CPZ™	1
EZND-100	100 mL EZ-CPZ-ND™	1
EZND-500	500 mL EZ-CPZ-ND™	1
STORZ-100	100 mL StorZ™	1
STORZ-500	500 mL StorZ™	1
STORZ-ND-100	100 mL StorZ-ND™	1
STORZ-ND-500	500 mL StorZ-ND™	1

EZ-CPZ™, EZ-CPZ-ND™, StorZ™, and StorZ-ND™ are trademarks of INCELL Corporation

All data and statements concerning these products may be considered as being indicative of representative properties and characteristics obtainable. We make no warranty, express or implied, concerning actual use or results because of industry specific influences.

