

CT-014-00

CT-014-90

CT-014-CVT

## Transfer with Confidence

- » For sterile transfer of cell based products, including whole blood and cell suspensions, from one container to another
- » Designed to support cell therapy manufacturing
- » Compatible with DMSO up to 20%
- » Manufactured without DEHP
- » Cleanroom compatible packaging
- » Compatible with sterile connection technology

## Ordering Information

### Standard Applications

- » Sterile transfer of whole blood and cell suspensions

Catalog #	Description	Case Qty.
CT-014-00	24" PVC with 2 spike adapters	60
CT-014-90	4" PVC with spike adapter and luer	60
CT-014-92	2.25" PVC with spike adapter and luer	60
CT-014-94	2.25" PVC with spike and male luer	60
CT-014-50	Spike with female luer and twist off	60
CT-014-06	6 lead harness with port adapter	24
CT-014-08	8 lead harness with port adapter	24
CT-014-10	10 lead harness with port adapter	24
CT-014-XL	5mL syringe with PVC lead and spike adapter	60
CT-014-SSS	60mL syringe with 6" PVC lead	36
CT-014-SSV	12mL syringe with 3 PVC leads and adapters	60
CT-014-VAS	6" PVC lead with syringe and vial adapter	60
CT-014-SPS	3 PVC leads with spike, port, syringe adapters	60
CT-014-RPS	3 PVC leads with sampling valve and spike	60
CT-014-CVT	10" C-Flex® to PVC Lead	60
CT-014-UPS	12 PVC leads with adapters	24



Pending 510(k) submission, not currently available in USA

All data and statements concerning these products may be considered as being indicative of representative properties and characteristics obtainable. We make no warranty, express or implied, concerning actual use or results because of industry specific influences.

