


# Plasma/Fluid Transfer Sets

for blood or blood components



- » For safe and efficient transfer of blood or blood components
- » Needleless design prevents accidental needle sticks
- » Some designs compatible with sterile connection devices
- » May be used for Quality Control Sampling
- » Attaches to a flexible, non-vented container to prepare, sample, or process blood components
- » Aids in compliance with national needle stick safety laws
- » Alternative to existing sampling site couplers and potentially hazardous needle/syringe systems

## Standard Applications

- » Perform sterile transfer of blood or blood components from one container to another

## Ordering Information

Catalog #	Description	Case Qty.
03-220-00	24" tubing, roller clamp, and two piercing pins	96
03-220-30	24" tubing, roller clamp, and needle adapter	96
03-220-50	3" tubing with piercing pin, dual syringe adapter, and spike port adapter	96
03-220-90	4" tubing with piercing pin and syringe adapter (female luer)	96
03-220-92	2" tubing with piercing pin and syringe adapter (female luer)	96
03-220-94	2" tubing with piercing pin, pinch clamp, and male luer adapter	96
03-220-XL	6" tubing with piercing pin and removeable 5mL BD™ Syringe	96



All data and statements concerning these products may be considered as being indicative of representative properties and characteristics obtainable. We make no warranty, express or implied, concerning actual use or results because of industry specific influences.

