




# Blood Management

 charterMedical


## Single-use solutions for **blood** and **blood components**:

- Administration
- Collection
- Processing
- Separation
- Storage
- Transport



Visit [www.chartermedical.com/bloodmgmt](http://www.chartermedical.com/bloodmgmt)

Charter Medical Blood Transfusion Products are sold exclusively through a specialty distribution network to hospitals, plasma centers, blood centers, and healthcare providers in the US only.



03-960-00F

03-960-33

N40-960-33

## Neonatal Pediatric Syringe Sets

Catalog #	Description	Case Qty
03-960-00	Blood administration set w/ 150µm filter & 30cc BD™ syringe	24
03-960-00F	Blood administration set w/ 150µm filter & 30cc BD™ syringe - not manufactured w/ DEHP	24
03-960-10	Blood administration set w/ 150µm filter & 50cc BD™ syringe	24
03-960-10F	Blood administration set w/ 150µm filter & 50cc BD™ syringe - not manufactured w/ DEHP	24
03-960-10FX	Spikeless blood administration set w/ 150µm filter & 50cc BD™ syringe - not manufactured w/ DEHP	24
03-960-10X3	Blood administration set w/ 150µm filter & three 50cc BD™ syringes	24
03-960-32	Blood administration set w/ 150µm filter, side clamp & female luer lock adapter	48
03-960-33	Blood administration set w/ 150µm filter, pinch clamp & female luer lock adapter	48
03-960-33F	Blood administration set w/ 150µm filter, pinch clamp & female luer lock adapter - not manufactured w/ DEHP	48
03-960-34	Blood administration set w/ 150µm filter & spike port	48
03-960-35	Blood administration set w/ 150µm filter & 35cc MONOJECT® syringe	24
03-960-35F	Blood administration set w/ 150µm filter & 35cc MONOJECT® syringe - not manufactured w/ DEHP	24
03-960-37	Blood administration set w/ 150µm filter & 60cc MONOJECT® syringe	24
03-960-37F	Blood administration set w/ 150µm filter & 60cc MONOJECT® syringe - not manufactured w/ DEHP	24
03-960-37FX	Spikeless blood administration set w/ 150µm filter & 60cc MONOJECT® syringe - not manufactured w/ DEHP	24
03-960-62	Blood administration set w/ 150µm filter & luer adapters	48
03-BPD-30	Blood administration set w/ 30cc BD™ syringe, attached 150mL bag & 150µm filter	18
03-BPD-60	Blood administration set w/ 50cc BD™ syringe, attached 150mL bag & 150µm filter	18
03-MJP-60	Blood administration set w/ 60cc MONOJECT® syringe, attached 150mL bag & 150µm filter	18
N40-960-00	NEO-40 blood administration set w/ 40µm filter & 30cc BD™ syringe	24
N40-960-10	NEO-40 blood administration set w/ 40µm filter & 50cc BD™ syringe	24
N40-960-33	NEO-40 blood administration set w/ 40µm filter, pinch clamp & female luer lock adapter	48
N40-960-35	NEO-40 blood administration set w/ 40µm filter & 35cc MONOJECT® syringe	24
N40-960-37	NEO-40 blood administration set w/ 40µm filter & 60cc MONOJECT® syringe	24



## Aliquot Bag Systems

Catalog #	Description	Case Qty
T2008	150mL eight-bag aliquot system w/ piercing pin (unsegmented tubing)	24
T2182	150mL three-bag aliquot system w/ piercing pin	24
T2186	150mL six-bag aliquot system w/ piercing pin	24
T2286	150mL six-bag aliquot system w/ piercing pin & needleless access port	24
T3000	150mL four-bag aliquot system w/ piercing pin	24
T3002	150mL two-bag aliquot system w/ piercing pin	24
T3003	150 mL four-bag aliquot system w/ 150µm filter piercing pin & needleless access port	24



03-220-50



03-220-90

## Plasma/Fluid Transfer Sets

Catalog #	Description	Case Qty
03-220-00	24" tubing w/ roller clamp & two piercing pins	96
03-220-30	24" tubing w/ roller clamp, piercing pin & slip lock (needle) adapter	96
03-220-50	3.5" inlet tubing w/ pinch clamp, piercing pin, dual syringe adapter & spike port adapter	96
03-220-90	4" tubing w/ pinch clamp, piercing pin & female luer lock (syringe) adapter	96
03-220-92	2" tubing w/ pinch clamp, piercing pin & female luer lock (syringe) adapter	96
03-220-94	2" tubing w/ pinch clamp, piercing pin & male luer lock adapter	96
03-220-XL	6" tubing w/ piercing pin & removable 5mL BD™ syringe	96




# Plasma/Platelet Sampling Devices

Catalog #	Description	Case Qty
03-220-100	8" Platelet sampling set with bacterial detection viral adapter cap	48
03-220-BC	Platelet sampling set with 10 mL BD™ syringe and cap/viral adapter	24
03-220-BC12	Platelet sampling set with 12 mL MONOJECT® syringe and cap/viral adapter	96
03-220-BC-20	Platelet sampling set with 20cc BD™ syringe and cap/viral adapter	20
03-220-TM	Platelet sampling set with 100 mL bag and safety needle assembly	48
03-220-TMBC	Platelet sampling set with 100 mL bag, 12 mL MONOJECT® syringe and cap/viral adapter	24



# Blood Transfer Bags

Catalog #	Description	Case Qty
T3101	150mL transfer bag w/ piercing pin	84
T3104	400mL transfer bag w/ piercing pin	48
T3106	600mL transfer bag w/ piercing pin	48
T3106X2	Two 600mL transfer bags w/ piercing pin	24
T3107	1000mL transfer bag w/ piercing pin	32
T3108	2000mL transfer bag w/ piercing pin	24
T3109	3000mL transfer bag w/ piercing pin	24
T3110	3000mL transfer bag w/ multiple luer adapters	40
C600	600mL autologous blood transfer pack w/ piercing pin	48
C605	600mL autologous blood transfer pack w/ needle adapter	48



T8006



T4105

## Pooling Harness Sets

Catalog #	Description	Case Qty
T3006	3000 mL bag with six-lead harness set, attached	24
T4004	400 mL bag with four-lead harness set, attached	24
T4005	400 mL bag with five-lead harness, attached	24
T4006	400 mL bag with six-lead harness, attached	24
T4105	400 mL bag with five-lead sealed tubing harness set, attached	24
T6010	600 mL bag with ten-lead harness, attached	36
T8006	Six-lead harness set with port adapter	36
T8010	Ten-lead harness set with port adapter	36



T3401

## Apheresis Plasma Processing Bags

Catalog #	Description	Case Qty
T3400	400mL triple bag aliquot system w/ slip lock adapter	24
T3401	300mL triple bag aliquot system w/ sealed tubing	24

# EXP-Pak™ Cell Expansion Containers

Catalog #	Description	Case Qty
EXP-50	30mL fill volume at 1 cm thick cell expansion bag w/ two leads & female luer, & dockable fill line w/ sampling valve	24
EXP-500	190mL fill volume at 1cm thick cell expansion bag w/ two leads & female luer, & dockable fill line w/ sampling valve	12
EXP-1L	380mL fill volume at 1cm thick cell expansion bag w/ two leads & female luer, & dockable fill line w/ sampling valve	12
EXP-3L	630mL fill volume at 1cm thick cell expansion bag w/ two leads & female luer, & dockable fill line w/ sampling valve	12
EXP-5L	1250mL fill volume at 1cm thick cell expansion bag w/ two leads & female luer, w/ sampling valve	12



Looking for a DEHP-free cell expansion option?  
Ask about our CXP containers.

# Cell Freeze® Cryogenic Storage Containers

Catalog #	Description	Case Qty
CF-50	10-20mL fill volume cryogenic storage container w/ two PVC leads w/ female luer	40
CF-250	30-70mL fill volume cryogenic storage container w/ two PVC leads w/ female luer	40
CF-250-3L	30-70mL fill volume cryogenic storage container w/ two PVC leads w/ female luer & one male luer	40
CF-500	55-100mL fill volume cryogenic storage container w/ two PVC leads w/ female luer	40
CF-500-3L	55-100mL fill volume cryogenic storage container w/ two PVC leads w/ female luer & one male luer	40
CF-750	80-190mL fill volume cryogenic storage container w/ two PVC leads w/ female luer	40



Contact the distributor in your region to order these critical healthcare solutions. Click [here](#) for a map of distributors to find the one that serves your area, or visit [www.chartermedical.com/bloodmgmt](http://www.chartermedical.com/bloodmgmt).

BD is a trademark of Becton, Dickinson and Company.  
MONOJECT is a registered trademark of Covidien AG or an affiliate.  
Cell Freeze and EXP-Pak are trademarks and/or registered trademarks of Charter Medical, LLC.

All data and statements concerning these products may be considered as being indicative of representative properties and characteristics obtainable. We make no warranty, express or implied, concerning actual use or results because of industry specific influences.

