


# Apheresis Plasma Processing Bags

## for Plasma Donation

- For aliquoting a jumbo plasma unit into multiple plasma units
- Regulated medical device

### Benefits

- Compatible with sterile connection technology
- Compatible with commercial plasma apheresis technology to facilitate plasma collection
- Available in multiple configurations



T3401

**Intended Use Statement**

For collection, storage, processing, transport, separation, and administration of blood and blood components.

## Ordering Information

Catalog #	Description	Case Qty
T3400	400 mL triple bag aliquot system with slip lock adapter	24
T3401	300 mL triple bag aliquot system with RF-welded tubing	24

Contact us today at [www.chartermedical.com](http://www.chartermedical.com)

All data and statements concerning these products may be considered as being indicative of representative properties and characteristics obtainable. We make no warranty, express or implied, concerning actual use or results because of industry-specific influences.