



Neonatal Syringe Set

with 40-micron Filter

- Needleless design to prevent accidental needle sticks, which complies with OSHA requirements and national needle stick safety laws¹
- Provides final filtration in accordance with Association for the Advancement of Blood & Biotherapies blood administration requirements²
- Regulated medical device

Benefits

- Provides closed blood-delivery system
- Maintains primary unit expiration date
- Maximizes usage of adult-size blood units
- Limits neonatal donor exposure
- Reduces aliquot preparation costs

Usage Details

- Syringe blood aliquot (RBC, Plasma) must be administered to patient within 24 hours and syringe platelet aliquot within 4 hours³
- Filter hold-up volume is less than 5 mL



Intended Use Statement

To remove microaggregates and particulate debris from whole blood and red blood cells during infusion

Configurations Available

- BD™ or MONOJECT® Syringes with 30cc, 35cc, 50cc, and 60cc syringe volumes

Ordering Information

Catalog #	Description	Case Qty
N40-960-00	Blood administration set with 40-micron filter and 30cc BD™ syringe	24
N40-960-10	Blood administration set with 40-micron filter and 50cc BD™ syringe	24
N40-960-33	Blood administration set with 40-micron filter and syringe adapter	48
N40-960-35	Blood administration set with 40-micron filter and 35cc MONOJECT® syringe	24
N40-960-37	Blood administration set with 40-micron filter and 60cc MONOJECT® syringe	24

Contact us today at www.chartermedical.com

References:

1. The Federal Needlestick Safety and Prevention Act (H.R. 5178) mandates needle stick safety practices nationwide.
2. Standards for Blood Banks and Transfusion Services. Bethesda, MD: AABB.
3. Standards for Blood Banks and Transfusion Services. Bethesda, MD: AABB.

BD is a trademark of Becton, Dickinson and Company.

MONOJECT is a registered trademark of Covidien AG or an affiliate.

All data and statements concerning these products may be considered as being indicative of representative properties and characteristics obtainable.

We make no warranty, express or implied, concerning actual use or results because of industry-specific influences.



N40SS_II/2024_VER3

Charter Medical is ISO 13485:2016 certified and FDA registered.

3948-A Westpoint Blvd., Winston-Salem, NC 27103, USA | 866.458.3116 | www.chartermedical.com

