

Blood Transfer Bags

For blood transfer and storage

- Designed for the collection, storage, processing, transport, separation, and administration of blood and blood components.

Benefits

- Ergonomically shaped spike designed for ease of use
- Easy grip tab for port access
- Tamper-evident spike port
- Printed segment tubing numbers for positive ID of retention segments
- Rounded shoulder design ensures maximum product recovery



Standard Applications

Designed for the collection, storage, transport, separation, and administration of blood and blood components

Other Blood Processing Solutions

Charter Medical offers additional solutions such as aliquot bag systems, blood transfer bags, cryogenic storage containers, neonatal syringe sets, platelet sampling sets, pooling harnesses, and plasma/fluid transfer solutions

Ordering Information

Catalog #	Description	Case Qty
T3101	150mL transfer bag with piercing pin	84
T3104	400mL transfer bag with piercing pin	48
T3106	600mL transfer bag with piercing pin	48
T3106X2	2- 600mL transfer bags with piercing pin	24
T3107	1000mL transfer bag with piercing pin	32
T3108	2000mL transfer bag with piercing pin	24
T3109	3000mL transfer bag with piercing pin	24
T3110	3000mL transfer bag with 2 male and 1 female luer connector	40
C600	600mL autologous transfer bag with piercing pin	48
C605	600mL autologous transfer bag with needle adapter	48

Contact us today at www.chartermedical.com/contact-us

All data and statements concerning these products may be considered as being indicative of representative properties and characteristics obtainable. We make no warranty, express or implied, concerning actual use or results because of industry-specific influences.