


Platelet Sampling Devices

For safe and efficient transfer of platelet samples

- Designed for performing a sterile transfer of blood or blood components from one container to another.

Benefits

- Custom tubing compatible with conventional sterile docking devices
- Easy-to-use designs ensure sterile sampling of platelets for bacterial testing
- Designed to deliver precise volume of platelets to culture bottles
- Integrated bioMérieux culture vial adapter cap



Standard Applications

Designed for performing a sterile transfer of blood or blood components from one container to another

Other Blood Processing Solutions

Charter Medical offers additional solutions such as aliquot bag systems, blood transfer bags, cryogenic storage containers, neonatal syringe sets, pooling harnesses, and plasma/fluid transfer solutions

Ordering Information

Catalog #	Description	Case Qty
03-220-BC	10mL BD™ syringe and cap/vial adapter	24
03-220-BC12	12mL syringe and cap/vial adapter	96
03-220-TM	100mL bag and safety needle assembly	48
03-220-TMBC	100mL bag, 12mL syringe, and cap/vial adapter	24
03-220-BC-20	20mL BD™ syringe and cap/vial adapter	20

Contact us today at www.chartermedical.com/contact-us

All data and statements concerning these products may be considered as being indicative of representative properties and characteristics obtainable. We make no warranty, express or implied, concerning actual use or results because of industry-specific influences.